## Michigan Collaborative Standardization of Compounded Oral Liquids

Drug name:	Disopyramide	Dosage Form:	Suspension
<b>Concentration:</b>	10 mg/mL	Shelf Life:	28 days
Route:	Oral	Storage:	Room Temp
Volume:	100 mL	Auxiliary Labeling:	Shake Well

Ingredients	QS	Quantity	Units
Disopyramide 100 mg		10	Capsules
Cherry syrup	Y	100	mL

## **Directions:**

- 1. Open disopyramide capsules and empty contents into mortar
- 2. Triturate contents to a fine powder
- 3. Add a small amount of cherry syrup and levigate to a paste
- 4. Add cherry syrup via geometric dilution almost to volume, mixing thoroughly after each addition
- 5. Transfer contents of mortar to graduated cylinder
- 6. Rinse mortar and pestle with cherry syrup and pour into graduated cylinder
- 7. QS to 120 mL with cherry syrup
- 8. Transfer to glass amber bottle
- 9. Shake well

## **References:**

- Mathur LK et al. Am J Hosp Pharm. 1982;39:309-10.
- Jew RK, Soo-Hoo W, Erush SC. *Extemporaneous Formulations for Pediatric, Geriatric, and Special Needs Patients.* 2<sup>nd</sup> ed. Bethesda, MD: American Society of Health-System Pharmacists; 2010.

Reviewed: April 2013; ME/JR

The information provided on this site is only to be used for educational and informational purposes only. It may not be relied upon to direct or provide care to a patient, or to prepare treatments for a patient. This material is not intended as a substitute for retaining and consulting appropriate medical and health care professionals, who remain solely responsible to provide and direct patient care. Be aware that the information contained on this site may not be up to date due to more recent medical developments. The preparation, dispensing, and administration of pharmaceuticals is and remains the responsibility of licensed medical professionals or pharmacists, and may require governmental approvals. The owner and provider of this site disclaims and assumes no liability for any injury and/or damage to persons from use of any information contained on this site.

**ERSITY OF MICHIGAN** 





